INSTRUCTIONS FOR USE

KTR Nichro Stainless Steel Primary Molar Crowns

Instructions for Use, Sterilization Guidance, and MRI Compatibility Information

1. Product Description

KTR Nichro Stainless Steel Primary Molar Crowns are specifically crafted to provide durable provisional protection for pediatric posterior dentition. These crowns are pre-sized, pre-contoured, and pre-cripped to reduce the need for significant trimming and reshaping during placement. Additionally, they incorporate a progressively thinner cervical wall to promote optimal marginal fit.

Material

Crowns are manufactured from medical-grade 300-series (Type 304) nickel-chromium stainless steel. Size markings are designed to remain legible after normal sterilization cycles.

Intended Use

Durable provisional full-coverage protection of primary posterior teeth (primary molars) in pediatric patients. For use by licensed dental professionals only.

WARNINGS AND PRECAUTIONS

- Stainless steel alloy contains both nickel and chromium. A small portion of patients may exhibit hypersensitivity to these metals. Should any allergic signs present, advise the patient to consult a medical professional. Do not use in patients with a known hypersensitivity to nickel or chromium.
- Once a crown is permanently cemented it is a single-patient device and must never be reused.
- Any crown tried in a patient's mouth must be fully cleaned and sterilized before it is used on another patient (see Section 3).
- Stainless steel crowns can produce image artifacts on head and neck MRI scans (see Section 4).

2. Placement Instructions

Two procedures are described below. Procedure A applies when the selected crown requires no extensive custom adjustment. Procedure B applies when custom modification (trimming, recontouring, or crimping) is necessary before cementation.

A. Placement When No Extensive Custom Adjustment Is Needed

- Select the Appropriate Crown Size:** Measure the mesio-distal dimension intraorally using a Boley gauge. Once the correct measurement is determined, match the measurement with the crowns in the kit to ensure it replicates the natural tooth dimension.
- Verify Dimensions:** Cross-check with any preoperative measurements taken prior to preparation.
- Try-In:** Seat the selected crown on the prepared tooth and evaluate its occlusion relative to adjacent and opposing teeth. For mandibular teeth, start seating lingually and then apply pressure buccally. For maxillary teeth, seat the buccal side first, followed by the lingual. Use an explorer to inspect the margin and adjust if necessary to eliminate any small gaps.
- Reassess Occlusion:** Remove the crown and verify the occlusal contacts between the arches. If the crown interferes with proper occlusion, perform additional reduction and re-test the fit.
- Isolate and Prepare:** Rinse and isolate the prepared tooth using cotton rolls. Clean and dry the preparation thoroughly.
- Mix and Load Cement:** Prepare the luting agent as directed by the manufacturer and fill the crown, ensuring no air bubbles become trapped.
- Seat the Crown:** Follow the seating technique from Step 3. Allow excess cement to extrude naturally under pressure; avoid forceful final seating with excessive thumb pressure.
- Check Final Fit:** Quickly assess the crown's position relative to neighboring teeth and verify occlusion.
- Have Patient Bite to Seat:** Instruct the patient to close into normal occlusion to ensure accurate seating. When using a bite stick or tongue blade, distribute force across the quadrant mesio-distally. Never concentrate force only between the crown and its opposing tooth, as this may drive the crown too deep and damage the gingiva.

- Clean Excess Cement:** After the recommended setting time, remove any residual cement with an explorer or scaler. Use knotted floss to clean interproximal areas and rinse the gingival sulcus to remove debris.
- Evaluate Final Result:** The crown should closely replicate the natural tooth in arch space, contacts, and functional occlusion. A radiograph should confirm complete seating and a precise marginal fit without gaps.

B. Placement When Custom Modification Is Necessary

- Select the Best-Fitting Crown:** Choose a crown from the kit that matches the tooth's mesio-distal measurement (refer to Procedure A, Step 1).
- Confirm with Measurements:** Compare with any preoperative dimensions to ensure accuracy.
- Try-In:** If the crown causes excessive gingival blanching when seated, slight trimming may be necessary.
- Mark Trim Line:** Outline the gingival contour on both buccal and lingual surfaces to guide trimming.
- Trim to Length:** Use crown scissors to cut about 1 mm below the marked line. Blend edges to avoid sharp burrs that could irritate tissue; smooth with a heatless stone if needed.
- Recontour Margins:** Use contouring pliers to reshape any areas that were trimmed, ensuring a snug fit.
- Test Fit Again:** Seat the crown to confirm adaptation.
- Check Contacts:** Use floss to ensure proper contact points remain. Inspect margins with an explorer for any gaps.
- Verify Occlusion:** Ensure there are no occlusal interferences.
- Refine Fit for Snap Seating:** If contact points need reestablishing, adjust interproximal walls with contouring pliers.
- Final Crimping:** Adapt the margins all around the crown with crimping pliers to secure a tight seal.
- Smooth and Polish Edges:** Thin and polish the margins with a heatless stone to prevent soft-tissue irritation and to facilitate easier cleaning for the patient. Perform a final trial fit before cementation.
- Proceed with Cementation:** Follow Procedure A, Steps 5–11, to complete the final cementation process.

3. Cleaning and Sterilization Instructions

General Information

- KTR pediatric stainless steel crowns are made from medical-grade 300-series (Type 304) nickel-chromium stainless steel and can safely go through routine steam sterilization (autoclave).
- Size markings are designed to stay readable after normal sterilization cycles.
- Once a crown is cemented permanently, it is single-patient only and must never be reused.
- Any crown tried in a patient's mouth must be fully cleaned and sterilized before being used on another patient.
- Follow your autoclave manufacturer's validated cycles, local regulations, and CDC dental sterilization guidelines.
- Preferred method:** steam sterilization (autoclave) at 121–132 °C.

Key Principles

- Classification:** KTR crowns become "critical items" once placed in the mouth. Any tried-in crown must be sterilized before reuse.
- Do not reuse:** Crowns that have been permanently cemented.
- Preferred sterilization method:** Steam sterilization (autoclave).
- Avoid:** Metal brushes, steel wool, and abrasive pads. Glutaraldehyde or other cold sterilants.

Step-by-Step Sterilization Protocol

SCOPE

- This protocol applies to crowns that are **not** cemented and are being reprocessed for reuse (for example, crowns that were tried in but not placed).

A. Point-of-Use Handling

- Immediately after try-in:** Wipe off debris with water-moistened gauze. Place in a covered, leak-proof container for contaminated instruments.
- Prevent drying:** If cleaning will be delayed, keep crowns in water or an enzymatic holding solution.

B. Cleaning (Required Before Sterilization)

Cleaning is essential. Sterilization will not be effective if debris remains.

- 1. Pre-rinse:** Rinse under running water.
- 2. Detergent cleaning:** Use an ultrasonic cleaner with an enzymatic or neutral detergent, or a washer-disinfector approved for dental instruments. Ultrasonic tips: fully immerse the crowns, avoid piling, and follow the manufacturer's time (typically 5–10 minutes).
- 3. Avoid abrasives and harsh chemicals:** No metal brushes, steel wool, chlorine, or glutaraldehyde.
- 4. Rinse thoroughly:** Remove all detergent residues.
- 5. Dry:** Use lint-free towels or filtered air. Make sure crowns are fully dry before packaging.

C. Inspection and Sorting

- 1. Inspect under good lighting:** Look for debris, stains, cracks, perforations, corrosion, and legible size markings.
- 2. Discard damaged crowns:** Do not reuse crowns with defects.
- 3. Organize:** Sort by size and tooth type in labeled cassettes or pouches.

D. Packaging

- 1. Use FDA-cleared wraps or pouches:** Place crowns in cassettes or sterilization pouches.
- 2. Include indicators:** An external indicator on the pack and an internal indicator inside the pack.
- 3. Label:** Sterilizer ID (if multiple units), cycle/load number, date, and optional operator initials and crown size.

E. Steam Sterilization (Autoclave)

Use the cycle validated by your autoclave manufacturer. Typical CDC-referenced wrapped cycles are shown below.

Cycle type	Temperature	Exposure time
Gravity displacement	121 °C (250 °F)	30 minutes
Pre-vacuum	132 °C (270 °F)	5 minutes

- 1. Load correctly:** Do not overload. Place pouches on edge, paper side up.
- 2. Run the proper cycle:** Select the validated wrapped-instrument cycle.
- 3. Monitor:** Mechanical: time, temperature, pressure. Chemical: indicator color change. Biological: weekly spore testing per CDC.
- 4. Dry:** Allow full drying inside the autoclave.

F. Cooling, Storage, and Handling

- 1. Cooling:** Let packs cool to room temperature before handling.
- 2. Check packaging:** Look for tears, holes, wet spots, or failed indicators. Reprocess if compromised.
- 3. Storage:** Keep in a clean, dry, closed cabinet. Follow “event-related sterility” — packs stay sterile unless damaged or wet.
- 4. Opening at chairside:** Open aseptically and avoid touching the crown; use sterile forceps or drop the crown onto a sterile tray.

KTR Crowns is committed to supporting safe, consistent, and compliant infection-control practices. Following these guidelines helps ensure every patient receives the highest standard of care while protecting your team and preserving the quality of your stainless steel crowns.

4. MRI Compatibility Information

Pediatric stainless steel crowns (SSCs) are generally considered safe with respect to the physical risks of magnetic resonance imaging (MRI), but they are a known source of image interference. The primary concern is not that the crowns will dislodge or heat, but that they can create “blind spots” in the scan that might lead a physician to miss an important diagnosis.

Safety vs. Image Quality

Current clinical consensus and studies differentiate between the physical safety of the patient and the diagnostic quality of the image:

- Physical safety (low risk):** KTR SSCs are made of 300-series stainless steel (specifically Type 304). While this alloy contains nickel and chromium, it is largely non-ferromagnetic or only weakly magnetic. It is therefore unlikely to be pulled off the tooth or to heat significantly during a standard MRI.
- Image artifacts (high risk):** This is the main issue. Stainless steel creates “susceptibility artifacts” — distortions that appear as black holes or wavy lines on the scan.

- In-vivo results:** A 2024–2025 systematic review confirmed that SSCs are the only restorative material consistently associated with significant artifacts on head MRIs in children.
- Scope of interference:** These artifacts can obscure the jaw, tongue, and base of the skull. While they usually do not compromise scans of the brain's upper regions, they can interfere with views of the pituitary gland or the auditory nerves.

Key Study Findings

- Systematic reviews:** A recent review in the *European Archives of Paediatric Dentistry* noted that while there is no evidence of patient injury, the diagnostic value of a head and neck MRI can be compromised by the presence of these crowns.
- Material comparison:** Studies comparing SSCs to zirconia (ceramic) crowns found that zirconia is entirely MRI-compatible and produces almost no distortion. For children with known chronic conditions requiring frequent MRIs (for example, oncology or neurology patients), zirconia is now strongly recommended over stainless steel.
- Multiple crowns:** American Academy of Pediatric Dentistry (AAPD) and radiology guidance suggests that if a child has multiple SSCs (for example, four or more), the “magnetic load” is higher, increasing the size of the image distortion.

Summary of MRI Impacts

Concern	Risk level	Details
Movement / dislodgement	Very low	Modern dental cements are strong enough to withstand the weak magnetic pull of 304-grade steel.
Heating / burns	Very low	Small dental restorations rarely absorb enough radiofrequency energy to cause thermal injury.
Image distortion	High	Can create significant “shadows” on head, neck, and brain scans.

Guidance for Parents and Caregivers

Dental professionals may share the following guidance with parents and caregivers of patients who have stainless steel crowns:

- Notify the technologist:** Always tell the MRI technologist that the child has stainless steel crowns. They can often adjust the pulse sequences (the way the machine acquires the image) to minimize metal distortion.
- Location matters:** If the MRI is of a knee, spine, or abdomen, the dental crowns will have no impact on scan quality.
- Future planning:** If the child is likely to need repeated brain or head scans (for example, due to seizures or tumors), discuss with the pediatric dentist whether replacing SSCs with zirconia or composite crowns is appropriate.

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